



# DEPPERLER

*Signature of excellence*

# SMARTENDO KIT

The ultimate upgrade  
for essentials instruments  
in your daily work

## Prexo M

- Plug GP in apical region with control
- Apply biomaterials

## Explora

- Canal scouting at any level
- Manage intracanal obstacles


## MisurEndo

- Diagnose safety
- Smart Guide for the correct access cavity design

## ExtEndo\*

- Extend the length of GP cones in longer canals
- Remove GP from the canals

## Mirror



 Treatment & retreatment

 MB2 scouting

 Broken file removal

 MTA

*\*With heavy lateral pressure, the instrument is likely to break. We recommend replacing the guttapercha in the middle of the canal, and then scratching out along the wall.*

## Multiplo

- Measure files
- Adjust the GP cone length
- Guide GP cone and collect sealers



**DEPPERLER**

*Signature of excellence*

# WHY CHOOSE SMARTENDO?

Discover  
articles and  
clinical cases



*Efficiency*

**the swiss knife for 95% of daily endo work**

*Safety*

**with shapes and angulations being redesigned by endo masters for safe work**

*Ergonomy*

**thanks to steel quality, hand made process and instruments shapes**

*Clarity*

**with color coding for better instrument identification**



Developed in collaboration with  
Style Italiano Endodontics group.  
Talented specialists who share teachable,  
feasible and repeatable knowledge about  
endodontics.

[endodontics.styleitaliano.org](http://endodontics.styleitaliano.org)