




DEPPELER

Signature of excellence

GRACEY CURETTES

Area-specific curettes
for deep subgingival scaling



Gracey 1/2 | 1GE2CN



Gracey 3/4 | 3GE4CN



Gracey 5/6 | 5GE6CN



GRACEY CURETTES

THE CLASSIC RANGE

- High quality materials
- Very efficient cutting edges
- Excellent tactile feedback
- Ergonomic and comfortable

Gracey 7/8 | 7GE8CN



Gracey 9/10 | 9GE10CN



Gracey 11/12 | 11GE12CN



Gracey 13/14 | 13GE14CN



Gracey 15/16 | 15GE16CN



Gracey 17/18 | 17GE18CN



"No other system can rival these hand instruments in terms of precision, sturdiness and ease of handling."
Prof. Dr. Michael Haas,
Medical University of Graz, Austria

GRACEY CURETTES FOR ALL ANATOMIES

THE VARIATIONS

To complete the Classic range and meet specific needs, Deppeler created four variations of the 5/6, 7/8, 11/12 and 13/14 Gracey curettes: Deep, Mini, Nano and Rigid.



DEEP

ref. 5GL6CN, 7GL8CN, 11GL12CN, 13GL14CN
A longer shaft and adapted angulation for better access to deep periodontal pockets, while offering more comfort for the dentist and the patient.



MINI

ref. 5PM6CN, 7PM8CN, 11PM12CN, 13PM14CN
A shorter and thinner blade, optimal shaft and blade curvature to closely match root profile and irregularities, and minimize the risk of damaging soft tissues in narrow and deep pockets.



NANO


ref. 5HG6CN, 7HG8CN, 11HG12CN, 13HG14CN
A very short and thin blade, for gentle and easy access to all areas, designed to adapt to the smallest unevenness of the roots. Especially useful in the aftercare phase.



RIGID

ref. 5GC6CN, 7GC8CN, 11GC12CN, 13GC14CN
A larger shaft diameter to apply more pressure and remove heavy calculus on the tooth surface.

Identify your Gracey curettes at a glance with the handle color code.



Learn more
about the
Gracey curettes
range and
sharpening tips



Order information


These instruments can be ordered individually or as a kit.

Classic Gracey kit reference: **KITGECN**
includes

- ☒ Explorer – DH2
- ☒ Classic Gracey curette 5/6 – 5GE6CN
- ☒ Classic Gracey curette 7/8 – 7GE8CN
- ☒ Classic Gracey curette 11/12 – 11GE12CN
- ☒ Classic Gracey curette 13/14 – 13GE14CN
- ☒ Sterilization box – CLEANextB5

Handle options

CLEANext handle (fixed color code)



ref. 13GE14CN

ADEP grips (11 color options)



ref. 13GE14+ADEP B (color code)



WWW.DEPPERLER.CH

f DEPPELERSA | @ DEPPELER_INSTRUMENTS



DEPPERLER
Signature of excellence